

The World Leading Provider of High Pressure Equipment for Research and Industry since 1945!



Division of Snap-tite, Inc.

MAG075 Inline

MagneDrive® Series

At a Glance

Minimum Static Torque: 7 inch-lbs. (791 N-mm)

Material of


Construction: 316 Stainless Steel or Hastelloy® C276

Maximum Pressure: 6000 psi @ 850° F or 5000 psi @ 950° F
(414 bar @ 454° C or 344 bar @ 510° C)


Principle

The MagneDrive® agitator uses rare earth magnets, permitting packless mixing at higher speeds and with higher viscosity fluids. Outer drive magnets, rotated by a direct coupled motor exert powerful attraction on the encapsulated inner magnet assembly. As the outer drive magnets are rotated, the inner magnets are actuated, resulting in rotation of the agitator shaft.


The MagneDrive® Principle



External driver magnets



Encapsulated driver magnet assembly and sealed rotor shaft



Outer magnets are rotated by a motor, thus rotating inner magnets and rotor shaft.



Application and Benefits

The MagneDrive® Agitator is recognized worldwide as a highly efficient method of promoting chemical reactions and catalyst testing among gases, liquids and solids in high pressure autoclaves. It can be mated with any number of optional impellers, including our Dispersimax® turbine type gas dispersion impeller or with any spinning catalyst reactor baskets that are housed within a pressurized vessel. Custom engineering of impeller designs can be performed based on developed horsepower, viscosity, critical speed and other key factors associated with specific processes that need contamination-free, pack-less agitation.

Contamination-free mixing- Packless design eliminates shaft packing and need for lubrication.

Zero leakage to atmosphere- The MagneDrive® is a sealed system, closed to the atmosphere, so even sensitive fluids can be processed safely.

Continuous, high speed operation- No need to shut down in mid-reaction to change failed packing.

Liquid Cooled- Water cooling (user supplied) for over-temperature protection of magnets and bearings. Cooling flow is not always required and can vary depending on vessel operating temperature and drive speed.

Features

- Capable of mixing vessel sizes from 100 ml up to 4000 ml.
- Operating pressures as high as 6000 psi @ 850° F (414 bar @ 454° C).
- Direct in-line motor eliminates belts, reduces size, and creates nearly silent operation.
- Compact design with up to 7 in-lb (791 N-mm)³ of static torque.
- Designed for simple disassembly and maintenance. Bearings can be replaced with minimal effort.
- Carbon graphite (Purebon®) and Fluoropolymer with Carbon Fiber (FPCF) bearings available.
- Motors available up to 1/2 Hp².
- Various impellers available separately. Contact factory for details.

General Specifications

Maximum Speed: 3300 rpm¹

Material of Construction: All wetted parts 316 SS or Hastelloy C-276. For information on other materials, please consult factory.

Bearing Material: Purebon® 658RCH.⁴, Purebon® 3310.⁴ or Fluoropolymer with Carbon Fiber.

Maximum Pressure at Connection: 6000 psi @ 850° F or 5000 psi @ 950° F
(414 bar @ 454° C or 344 bar @ 510° C)

Maximum Temperature at Magnet Zone: 300° F (149°C)⁵

Maximum Temperature at Connection: 950° F (510°C)

Tachometer Pick-up: Magnetic hall effect sensor, pulse output

Shaft and Impeller: Supplied with standard length shafts. Customization is available. Autoclave Engineers offers a wide selection of impellers including the Dispersimax® gas dispersion system. Please consult factory for more information.

Liquid Cooling: Required for over-temperature protection of magnets and bearings. User supplied, standard water service, 0.3 gpm (1.1 lpm) min. flow rate.

Liquid Cooling Connections: 1/4" copper tube inlet and outlet.

¹ Maximum speeds may be limited by motor mixing requirements and shaft vibration, including critical speed.

² Motor horsepower should be sized at least 25% higher than the intended application requirement.

³ To determine horsepower at a certain speed, use the formula:

$$\text{hp} = \frac{T \times n}{63,025} \quad \text{where: } T = \text{torque in inch-lbs} \\ n = \text{speed in rpm}$$


⁴ Purebon is a registered Trademark of Pure Carbon.

⁵ The magnets are stabilized at 300° F (149° C). When the temperature of the magnets exceeds the stabilizing temperature for an extended period, loss of magnetic torque will occur. Some of this loss is not reversible and torque will not regenerate.

Please refer to the following sections of the catalog for complimentary products and additional technical details. See the *MAG075 Ordering Guide on the back cover to configure a drive for your specific application.*

| MAG075 Inline Magnedrive Drawings | Drawing Number |
|-----------------------------------|----------------|
| 316 Stainless Steel | 40C-0513 |
| Hastelloy® C-276 | 40C-0998 |
| Cover Connection | 10C-7227 |

| Magnedrive Dimesions | | | | | |
|---------------------------|----------------|------|--------------|-------|--------------------|
| Motor Type | Motor Dia. (Ø) | | Drive Height | | Shaft Length |
| | inches | mm | inches | mm | |
| 1/8 HP DC General Purpose | 3.41 | 86.6 | 15.09 | 383.3 | See Ordering Guide |
| 1/3 HP DC General Purpose | 3.41 | 86.6 | 17.32 | 439.9 | |
| 1/2 HP Air Motor | 3.65 | 92.7 | 13.82 | 351.0 | |



1/8 HP drive with Autoclave Engineers connection and typical shaft shown. Consult drawings (chart above) for additional dimensions.

Supporting Information

Dimensions

Ordering Guide

MAG075 A A B B C D E F G G G

| AA-Material | | | | | |
|--|---|-------------------------------------|-------|---------------------------|-------|
| SS | 316 SS | | | | |
| HC | Hastelloy ^{®1} C-276 | | | | |
| BB-Size | | | | | |
| 1I | 1" Magnet Stack (7 in-lb Static Torque) | | | | |
| C- Bearing | | | | | |
| 1 | Purebon [®] 658RCH ² | | | | |
| 2 | FPCF (Flouropolymer with Carbon Fiber) | | | | |
| 3 | Purebon [®] 3310 ² | | | | |
| D- Inline Motor | | | | | |
| D | 1/8 HP 0-130 VDC Variable Speed General Purpose (2500 rpm max.) | | | | |
| E | 1/3 HP 0-130 VDC Variable Speed General Purpose (2500 rpm max.) | | | | |
| F | Air Motor - Manual Speed Adjust (1/2 HP using 60 psi (4.1 bar) air @ 3000 rpm max.) | | | | |
| G | Air Motor - Electronic Speed Adjust 4-20 ma Input (1/2 HP using 60 psi (4.1 bar) air @ 3000 rpm max.) | | | | |
| E- Speed Sensor | | | | | |
| 0 | None | | | | |
| 2 | Intrinsically Safe Speed Sensor (No Barrier Provided) | | | | |
| F- Approvals | | | | | |
| 0 | None | | | | |
| C | CE Mark | | | | |
| GGG-Drive Shaft | | | | | |
| 000 | None (Shaft Purchased Separately) | | | | |
| See Table Below | Shaft Options for Standard Stirred Reactors | | | | |
| Vessel Style | Volume (ml) | Shaft Length (See Dimensional Fig.) | | Shaft Style Ordering Code | |
| | | Inches | mm | Dispersamax | Solid |
| Bolted Closure | 100 | 4.76 | 120.9 | D01 | S01 |
| | 300 | 8.51 | 216.2 | D02 | S02 |
| | 500 | 7.03 | 178.6 | D03 | S03 |
| | 1000 | 11.16 | 283.5 | D04 | S04 |
| EZE-Seal | 100 | 4.38 | 111.3 | D05 | S05 |
| | 300 | 8.09 | 205.5 | D06 | S06 |
| | 500 | 6.13 | 155.7 | D07 | S07 |
| | 1000 | 10.25 | 260.4 | D08 | S08 |
| Zipper Closure | 500 | 7.24 | 183.9 | D09 | S09 |
| | 1000 | 11.40 | 289.6 | D10 | S10 |
| Bolted Closure, EZE-Seal, Zipper Closure | 2000 | 8.63 | 219.2 | D11 | S11 |
| | 4000 | 14.88 | 378.0 | D12 | S12 |

NOTES:

¹HASTELLO[®] is a registered trademark of Haynes International Inc.

²Purebon[®] is a registered trademark of Morgan AM & T Inc.



Division of Snap-tite, Inc.

8325 Hessinger Drive
 Erie, Pennsylvania 16509-4679
 PH: 814-860-5700 FAX: 814-860-5811
 e-mail: aeclave_sales@snap-tite.com
www.autoclaveengineers.com

ISO-9001 Certified



06-1165SE-0611

! WARNING !

FAILURE OR IMPROPER SELECTION OR IMPROPER USE OF THE PRODUCTS AND/OR SYSTEMS DESCRIBED HEREIN OR RELATED ITEMS CAN CAUSE DEATH, PERSONAL INJURY AND/OR PROPERTY DAMAGE.
 This document and other information from Snap-tite, Inc., its subsidiaries and authorized distributors, provides product and/or system options for further investigation by users having technical expertise. It is important that you analyze all aspects of your application and review the information concerning the product or system in the current product catalog. Due to the variety of operation conditions and applications for these products or systems, the user, through its own analysis and testing, is solely responsible for making the final selection of the products and systems and assuring that all performance, safety and warning requirements of the application are met. The products described herein, including without limitation, product features, specifications, designs, availability and pricing, are subject to change by Snap-tite, Inc. and its subsidiaries at any time without notice.

Bulletin AGT-MAG075 Inline

The World Leading Provider of High Pressure Equipment for Research and Industry since 1945!